

SIRONIX

US Food and Drug Administration
FDA
510(k) Approved

STATIV[®]

US PATENTED ALL SUTURE SOFT ANCHOR IS NOW
USFDA APPROVED



MAY THE SPORT OF LIFE NEVER STOP


CONSTRUCTED FULLY OF STRONG ULTRA HIGH MOLECULAR WEIGHT POLYETHYLENE (UHMWPE) SUTURE MATERIAL,LOADED WITH SINGLE/DOUBLE/TRIPLE LOADED SUTURES(S)/TAPE(S).



A soft anchor with patented design provides benefits of bone preservation while giving fixation benefits of traditional anchors.



Anchor is loaded on minimalistic designed inserter device with bead to give assured centralized deployment of the anchor into pilot hole and in those difficult areas during labral repair & allows multiple anchor placement in glenoid rim.




The deployed unique tri-pod bunched up pattern doubles compared up to the original size of anchor hence gives much higher pull out strength*.




SHOULDER INSTABILITY



- Sleeves available for 1.5mm & 1.8mm anchors.
- Curved sleeve available with 12.5° and 25° angles.
- Drill bit with colour demarcation for 1.5mm & 1.8mm in black & red respectively.



ROTATOR CUFF REPAIR



- Strong dispenser.
- No sleeve or drill required.
- Available sizes - 2.5mm & 2.9mm.
- After insertion and punching, the tripod structure with more than 5.5mm.
- Recommend in Rotator cuff repair and Remplissage.

SHOULDER IMPLANTS ALL SUTURE ANCHOR

Catalogue No.	Description
S33-1501-SN	Stativ All Suture Anchor 1.5mm Neo-Single Loaded Suture
S33-15E1-SN	Stativ All Suture Anchor 1.5mm Neo-Single Loaded Tape
S33-1802-SN	Stativ All Suture Anchor 1.8mm Neo- Double Loaded Suture
S33-18E2-SN	Stativ All Suture Anchor 1.8mm Neo-Double Loaded Tape
S33-2502-ST	Stativ All Suture Anchor 2.5mm Self Tap Double Loaded Suture
S33-2503-ST	Stativ All Suture Anchor 2.5mm Self Tap Triple Loaded Suture
S33-25E2-ST	Stativ All Suture Anchor 2.5mm Self Tap Double Loaded Tape
S33-2902-ST	Stativ All Suture Anchor 2.9mm Self Tap Double Loaded Suture
S33-29E2-ST	Stativ All Suture Anchor 2.9mm Self Tap Double Loaded Tape
S33-2903-ST	Stativ All Suture Anchor 2.5mm Self Tap Triple Loaded Suture

SIRONIX ANCHOR INSTRUMENTATION

Catalogue No.	Description
S45-FL15-DU	1.5 mm Flexible Anchor Drill for Stativ 1.5mm curved cannula
S45-FL18-DU	1.8 mm Flexible Anchor Drill for Stativ 1.8mm curved cannula
S45-FL25-DU	2.5 mm Flexible Anchor Drill for Stativ 2.5mm curved cannula
S45-1512-LT	12.5° Curved Cannula for Stativ 1.5mm Anchor
S45-1525-LT	25° Curved Cannula for Stativ 1.5mm Anchor
S45-1812-LT	12.5° Curved Cannula for Stativ 1.8mm Anchor
S45-1825-LT	25° Curved Cannula for Stativ 1.8mm Anchor

Reference: *Data on file

Disclaimer:

Responsibility for administering products mentioned and interpretation of clinical data provided lies with the registered medical practitioner.
For more information about products, their application & precautions: please visit: www.healthiummedtech.com



To know more scan through your camera:

Healthium Medtech Ltd.

RMZ North Star, Cowrks 12th Floor, Adjacent to RMZ Galleria Mall, Yelahanka, Bangalore - 560 064, India
+(91)-(80)-4186 8198 | healthiummedtech.com

Healthium-Medtech

HealthiumMedtech

@HealthiumL

HealthiumMedtech

Sironix Medtech

@sironix_arthroscopysolutions

STATIV_B1_V2_14JUL21