


# MAY THE SPORT OF LIFE NEVER STOP



**SIRONIX**  
*Arthroscopy Solutions*


**FEMORAL FIXATION IMPLANTS**






SIRONIX PROLOOP™ Adjustable Loop Button Fixation device helps ensure accurate fixation in ligament reconstructive surgery.







**Secure double locking Loop design\***

**Ultra-high-molecular-weight polyethylene (UHMWPE) suture-** provides good fixation and tensile strength

**Interwoven loop-** Proprietary weaving process ensures a good mix of strength and flexibility

**Flexibility-** The adjustable loop meets the surgeon's flexibility requirement for fixation

**Button-** Titanium alloy, meets ASTM standard F136





## VARIABLE



The device is used for fixation of ligaments during orthopedic reconstruction procedures such as anterior cruciate ligament (ACL), posterior cruciate ligament (PCL), medial collateral ligament & other ligament to bone fixation surgeries

## FIXED



## ORDERING INFORMATION

### 12mm INFILOOP® FIXED LOOP BUTTON

Product Code	Description
S24-1215-C	15mm ACL Button
S24-1220-C	20mm ACL Button
S24-1225-C	25mm ACL Button
S24-1230-C	30mm ACL Button
S24-1235-C	35mm ACL Button

### 14mm INFILOOP® FIXED LOOP BUTTON

Product Code	Description
S24-1435-C	35mm ACL Button
S24-1430-C	30mm ACL Button
S24-1425-C	25mm ACL Button
S24-1420-C	20mm ACL Button
S24-1415-C	15mm ACL Button

## PROLOOP™ ADJUSTABLE LOOP BUTTON

S24-1436-C	60mm Adjustable Loop ACL Button
------------	---------------------------------

### PROLOOP™ ULTRA

14mm button with  
60mm / 90mm  
adjustable loop

### PROLOOP™ MINI

8mm button  
with 100mm  
adjustable loop

### PROLOOP™ XL

20mm button  
with 90mm  
adjustable loop

# SIRONIX

*Arthroscopy Solutions*

Reference: \*Data on file

**Disclaimer:**


Responsibility for administering products mentioned and interpretation of clinical data provided lies with the registered medical practitioner.  
For more information about products, their application & precautions: please visit: [www.healthiummedtech.com](http://www.healthiummedtech.com)


**Healthium Medtech Pvt Ltd.**

RMZ North Star, Cowrks 12<sup>th</sup> Floor, Adjacent to RMZ Galleria Mall, Yelahanka, Bangalore - 560 064, India  
+(91)-(80)-4186 8198 | [healthiummedtech.com](http://healthiummedtech.com)

 [Healthium-Medtech](https://www.facebook.com/Healthium-Medtech)

 [HealthiumMedtech](https://www.linkedin.com/company/HealthiumMedtech)

 [@HealthiumL](https://twitter.com/HealthiumL)

 [HealthiumMedtech](https://www.youtube.com/HealthiumMedtech)