

SIRONIX
Arthroscopy Solutions

REGAIN POWER
WITH SHOULDER IMPLANT
RANGE*



MAY THE SPORT OF LIFE NEVER STOP

SHOULDER IMPLANTS

BANKART & ROTATOR CUFF REPAIR SOLUTIONS

STATIV® ALL SUTURE ANCHOR

- A soft anchor with patented design provides benefits of bone preservation while giving fixation benefits of traditional anchors.
- Constructed of strong Ultra High Molecular Weight Polyethylene (UHMWPE) suture material.
- This soft anchor is loaded with sutures(s)/tape(s) pre-assembled and loaded on a minimalistic insertion device with bead to give assured centralized deployment of the anchor into pilot hole.
- The deployed unique tri-pod bunched up pattern doubles up compared to the original size of the anchor hence gives much higher pull out strength.



Available Options

Size	1.5mm	1.8mm	2.5mm
Suture Loading	Single Suture	Double Suture	Double/Triple Suture
Tape Loading	Single Tape	Double Tape	Double Tape

VIPILOK® KNOTLESS ANCHOR

- A versatile knotless anchor comes preloaded on an inserter assembly.
- Intended for fixation of suture to bone/tissue fixation.
- Capable of taking multiple (4-6) sutures/tapes from the tip eyelets*.
- Anchors available in wedge & screw design.



Available Options

Size	2.4mm	4.8mm	5.5mm
Tread Type	Wedge	Screw & Wedge	Screw & Wedge
Anchor Material Type	βTCP & PEEK	βTCP & PEEK	βTCP & PEEK
Tip Material Type	PEEK & Titanium	Titanium	Titanium

SHOULDER IMPLANTS

BANKART & ROTATOR CUFF REPAIR SOLUTIONS

CEPTRE® KNOTTED ANCHOR

- A synergistic combination of unique design, variety of material choices and high strength of UHMWPE suture.
- The unique pattern of threads allows easy insertion and provides maximum fixation strength.
- The anchors are double or triple loaded with suture(s)/tape(s) to provide the best possible combination for superior repair strength.
- The titanium anchor with fully threaded design suitable especially for weak/osteoporotic bone.



Available Options

Size	4.8mm	5.5mm
Loading Type	Double Suture & Double Tape	Double Suture, Double Tape & Triple Suture
Material Type	VIPLK, PEEK & Titanium	VIPLK, PEEK & Titanium

SPYKE™ KNOTTED ANCHOR

- A shoulder anchor range to offer easy insertion and optimum fixation for instability repair.
- Wedge design allows the cancellous bone to get proper fixation once inserted fully, preventing pullback of anchors .
- The anchors are single or double loaded with suture(s) to provide the best possible combination for superior repair strength.
- The titanium anchor with fully threaded design suitable especially for weak/osteoporotic bone.



Available Options

Size	2.5mm	2.8mm
Loading Type	Single Suture	Double Suture
Material Type	VIPLK, PEEK & Titanium	VIPLK, PEEK & Titanium

SHOULDER INSTRUMENTS

S36-0001-T	Tissue grasper with Ratchet Straight	S36-0024-T	Rasp down bend 20 degree
S36-0003-T	Crochet Hook Straight	S36-0025-T	Tissue Liberator Blade Down
S36-0004-T	Switching Stick	S36-1006-T	Tape knot pusher
S36-0005-T	Calibrated Probe	S36-2000-T	R-Suture Passer 17mm Working Length
S36-0006-T	Knot Pusher Full Loop	S39-0011/0012-T	Cannula Introducer 8mm/6mm
S36-0007-T	Sliding Suture Cutter Straight	S39-2001-T	SP Needle With Lateral Opening
S36-0008/0009-T	Clever Hook Right/Left	S35-0031-T	Shoulder Set Box- small
S36-0010/0011-T	Sixer Left/Right	S36-3006/3007-T	Disposable Lasso 45 Degree Left / 45 Degree Right
S36-0012-T	Penetrating Grasper Straight	S45-0015-D	1.5 mm Anchor Drill for Stativ 1.5mm anchor
S36-0013-T	Rasp up bend 20 degree	S45-0018-DU	1.8 mm Anchor Drill for Stativ 1.8mm anchor
S36-0014-T	Tissue Liberator Blade Up	S45-0024-DU	2.4 mm Anchor Drill for Viplok Knotless 2.4mm PEEK anchor
S36-0018-T	Tape Cutter	S45-0025-DU	2.5 mm Anchor Drill for Spyke 2.5mm anchors
S36-0019-T	Suture Manipulator	S45-0028-DU	2.8 mm Anchor Drill for Spyke 2.8mm anchors
S36-0020-T	Mallet	S45-0075-LU	Universal Guiding cannula
S36-0021/0022-T	Shoulder Instrument Box-Large/Premium	S45-1048-LU	Blunt Straight Awl-4.8 mm for Cepter & Viplok PEEK Tip
S36-0023-T	Combo Grasper	S45-1055-LU	Blunt Straight Awl-5.5 mm for Ceptre & Viplok PEEK Tip

Reference: *Data on file

Disclaimer:

Responsibility for administering products mentioned and interpretation of clinical data provided lies with the registered medical practitioner.

For more information about products, their application & precautions: please visit: www.healthiummedtech.com

Healthium Medtech Pvt Ltd.

RMZ North Star, Cowrks 12th Floor, Adjacent to RMZ Galleria Mall, Yelahanka, Bangalore - 560 064, India
 +(91)-(80)-4186 8198 | healthiummedtech.com